


Anti-Tnf-A Rabbit Monoclonal Antibody For Wb, Ip, If/Icc, Elisa - P01375

Item Number: CM0005650



Introduction

High-specificity rabbit monoclonal antibody targeting human TNF- α (SWISS: P01375), valid for WB, IP, IF/ICC, and ELISA. Ideal for cytokine, apoptosis, and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	TNF- α Rabbit mAb
Remarks/Alias	DIF; TNFA; TNFSF2; TNLG1F; TNF-alpha; TNF- α
Species	Human
GeneID (Human)	7124
GeneID	7124
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 77-233 of human TNF- α (NP_000585.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P01375
Protein Weight	26 kDa
Shipping	Ice bag