


Anti-Tnf-A Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P01375

Item Number: CM0016171



Introduction

Rabbit polyclonal antibody targeting human TNF- α , validated for WB, IF-P, IHC-P, and ELISA. Reacts with human, mouse, and rat. Useful for studying inflammation, apoptosis, and immune regulation.

[Learn More](#)

Parameter	Value
Product Name	TNF- α Rabbit pAb
Remarks/Alias	DIF; TNFA; TNFSF2; TNLG1F; TNF-alpha; TNF- α
Species	Human
GeneID (human)	7124
GeneID	7124
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01375
Protein Weight	26 kDa
Shipping	Ice bag