


Anti-Tnf-A Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P01375

Item Number: CM0017885



Introduction

Rabbit polyclonal antibody against human TNF- α for WB, IF/ICC, IHC-P, ELISA. Detects human and mouse TNF- α . Suitable for inflammation, apoptosis, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	TNF- α Rabbit pAb
Remarks/Alias	DIF; TNFA; TNFSF2; TNLG1F; TNF-alpha; TNF- α
Species	Human
GeneID (Human)	7124
GeneID	7124
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 134-233 of human TNF- α (NP_000585.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P01375
Protein Weight	26 kDa
Shipping	Ice bag