

Anti-Tnf-Alpha Monoclonal Antibody For Wb, If/Icc, Elisa - P06804

Item Number: CM0004854

Introduction

CamelBio's TNF-alpha Rabbit mAb (mouse reactive) validated for WB, IF/ICC, ELISA. Targets the 26 kDa TNF-alpha cytokine involved in inflammation, apoptosis, and immune regulation.

[Learn More](#)

Parameter	Value
Product Name	TNF-alpha Rabbit mAb
Remarks/Alias	DIF; Tnfa; TNF-a; TNFSF2; Tnlg1f; Tnfsf1a; TNFalpha; TNF-alpha
Species	Mouse
Gene ID (Human)	7124
Gene ID	21926
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 80-235 of mouse TNF- α (NP_038721.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse
SWISS	P06804
Protein Weight	26 kDa
Shipping	Ice bag