

Anti-Tnf Rii/Tnfrsf1B/Cd120B Polyclonal Antibody For Wb, Elisa - P25119

Item Number: CM0023952

Introduction

Rabbit polyclonal antibody against mouse TNF RII/TNFRSF1B/CD120b (TNFRSF1B), validated for WB and ELISA. Recognizes human and mouse targets. Ideal for studying TNF signaling and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	TNF RII/TNFRSF1B/CD120b Rabbit pAb
Remarks/Alias	p75; TNFBR; Tnfr2; CD120b; TNF-R2; TNFR80; TNFRII; Tnfr-1; TNF-R75; TNF-R-II; TNF-alphaR2; TNFalpha-R2; TNF RII/TNFRSF1B/CD120b
Species	Mouse
Gene ID (Human)	7133
Gene ID	21938
Immunogen	Recombinant protein corresponding to a sequence within amino acids 23-258 of mouse TNF RII/TNFRSF1B/CD120b (NP_035740.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P25119
Protein Weight	50kDa
Shipping	Ice bag