

Anti-Tnf α Rabbit Polyclonal Antibody For Wb And Elisa - P16599

Item Number: CM0023636

Introduction

High-quality rabbit polyclonal antibody targeting rat TNF α , validated for Western blot and ELISA. Detects TNF-alpha in human, mouse, and rat samples. Ideal for research on inflammation, immunology, and cell death.

[Learn More](#)

Parameter	Value
Product Name	TNF α Rabbit pAb
Remarks/Alias	Tnfa; RATTNF; TNF-alpha; TNF α
Species	Rat
GeneID (human)	7124
GeneID	24835
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 80-235 of rat TNF α (NP_036807.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P16599
Protein Weight	26kDa
Shipping	Ice bag