

Anti-Tnfaip1 Polyclonal Antibody For Wb And Elisa - Q13829

Item Number: CM0024667

Introduction

Rabbit polyclonal antibody against human TNFAIP1 (gene ID 7126), a BTB/POZ adapter protein involved in RhoA ubiquitination. Suited for Western blot (predicted band 36 kDa) and ELISA, facilitating studies of cytoskeletal regulation and apoptosis.

[Learn More](#)

Parameter	Value
Product Name	TNFAIP1 Rabbit pAb
Remarks/Alias	B12; B61; EDP1; BTBD34; hBACURD2; TNFAIP1
Species	Human
GeneID (Human)	7126
GeneID	7126
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 90-160 of human TNFAIP1 (NP_066960.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q13829
Protein Weight	36kDa
Shipping	Ice bag