

Anti-Tnfaip3 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P21580

Item Number: CM0002399



Introduction

Validated anti-TNFAIP3 rabbit monoclonal antibody for Western blot, IHC paraffin, and ELISA. Detects human TNFAIP3/A20, a ubiquitin-editing enzyme critical in NF- κ B signaling and inflammatory responses. Ideal for immunology and cancer research.

[Learn More](#)

Parameter	Value
Product Name	TNFAIP3 Rabbit mAb
Remarks/Alias	A20; AISBL; AIFBL1; OTUD7C; TNFA1P2; TNFAIP3
Species	Human
Gene ID (Human)	7128
Gene ID	7128
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 450-550 of human TNFAIP3 (P21580).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P21580
Protein Weight	90kDa
Shipping	Ice bag