


Anti-Tnfaip3 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P21580

Item Number: CM0024052



Introduction

TNFAIP3 rabbit polyclonal antibody validated for WB, IHC-P, ELISA. Reacts with human and rat. Immunogen corresponds to a.a. 600-738 of human TNFAIP3. Targets the ubiquitin-editing enzyme A20, involved in NF- κ B signaling and inflammation. For research use.

[Learn More](#)

Parameter	Value
Product Name	TNFAIP3 Rabbit pAb
Remarks/Alias	A20; AISBL; AIFBL1; OTUD7C; TNFA1P2; TNFAIP3
Species	Human
GeneID (human)	7128
GeneID	7128
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 600-738 of human TNFAIP3 (NP_006281.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P21580
Protein Weight	90 kDa
Shipping	Ice bag