


Anti-Tnfaip6 Rabbit Polyclonal Antibody For Wb, Elisa - P98066

Item Number: CM0024908



Introduction

Rabbit polyclonal antibody against human TNFAIP6 (TSG-6). Suitable for WB and ELISA. Cross-reacts with human, mouse, and rat. Ideal for extracellular matrix, inflammation, and bone remodeling research.

[Learn More](#)

Parameter	Value
Product Name	TNFAIP6 Rabbit pAb
Remarks/Alias	TSG6; TSG-6; TNFAIP6
Species	Human
GeneID (Human)	7130
GeneID	7130
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-277 of human TNFAIP6 (NP_009046.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P98066
Protein Weight	31kDa
Shipping	Ice bag