


Anti-Tnfaip8L2 (Tipe2) Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - Q6P589

Item Number: CM0017960



Introduction

Anti-TNFAIP8L2 (TIPE2) Rabbit Polyclonal Antibody for human, mouse, rat WB, IF/ICC, IP, ELISA. TNFAIP8L2 negatively regulates innate/adaptive immunity; inhibits TCR-mediated T-cell activation and autolysosome formation. Ideal for immunology research.

[Learn More](#)

Parameter	Value
Product Name	TNFAIP8L2 Rabbit pAb
Remarks/Alias	TIPE2; TNFAIP8L2
Species	Human
GeneID (Human)	79626
GeneID	79626
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-184 of human TNFAIP8L2 (NP_078851.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6P589
Protein Weight	21kDa
Shipping	Ice bag