


Anti-Tnfr1/Tnfrsf1A Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P19438

Item Number: CM0013310



Introduction

High-quality rabbit polyclonal TNFR1/TNFRSF1A antibody for Western blot, IHC, IF/ICC, and ELISA. Detects human, mouse, and rat TNFR1. Essential tool for TNF receptor signaling and apoptosis studies.

[Learn More](#)

Parameter	Value
Product Name	TNFR1/TNFRSF1A Rabbit pAb
Remarks/Alias	FPF; p55; p60; TBP1; TNF-R; TNFAR; TNFR1; p55-R; CD120a; TNFR55; TNFR60; TNF-R-I; TNF-R55; TNFR1/TNFRSF1A
Species	Human
GeneID (Human)	7132
GeneID	7132
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 345-455 of human TNFR1/TNFRSF1A (NP_001056.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P19438
Protein Weight	50kDa
Shipping	Ice bag