


## Anti-Tnfrsf10A Polyclonal Antibody For Wb, Elisa - 000220

Item Number: CM0029628



### Introduction

Rabbit polyclonal antibody against human TNFRSF10A (DR4/TRAIL-R1). Recombinant immunogen (aa 269-468). Suitable for WB and ELISA applications. For apoptosis and death receptor signaling research.

[Learn More](#)

| Parameter        | Specification                                                                                                          |
|------------------|------------------------------------------------------------------------------------------------------------------------|
| Product Name     | TNFRSF10A Rabbit pAb                                                                                                   |
| Remarks/Alias    | DR4, APO2, CD261, TRAILR1, TRAILR-1, TNFRSF10A                                                                         |
| Species          | Human                                                                                                                  |
| GeneID (Human)   | 8797                                                                                                                   |
| GeneID           | 8797                                                                                                                   |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 269-468 of human TNFRSF10A (NP_003835.3) |
| Source           | Rabbit                                                                                                                 |
| Category         | Polyclonal Antibodies                                                                                                  |
| Application      | WB, ELISA                                                                                                              |
| Cross Reactivity | Human                                                                                                                  |
| SWISS            | O00220                                                                                                                 |
| Protein Weight   | 50kDa                                                                                                                  |
| Shipping         | Ice bag                                                                                                                |