


Anti-Tnfrsf10A Rabbit Polyclonal Antibody For Wb, If/Icc And Elisa - O00220

Item Number: CM0016096



Introduction

Rabbit polyclonal antibody targeting human TNFRSF10A (Death Receptor 4/TRAIL-R1). Confirmed for Western blot, immunofluorescence, and ELISA. Cross-reactivity with mouse and rat. Detects a band at approximately 50 kDa (glycosylated forms may appear larger). Ideal for investigating extrinsic apoptosis pathways.

[Learn More](#)

Parameter	Value
Product Name	TNFRSF10A Rabbit pAb
Remarks/Alias	DR4; APO2; CD261; TRAILR1; TRAILR-1; TNFRSF10A
Species	Human
Gene ID (Human)	8797
Gene ID	8797
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 269-468 of human TNFRSF10A (NP_003835.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00220
Protein Weight	50kDa
Shipping	Ice bag