

Anti-Tnfrsf10C Rabbit Polyclonal Antibody For Wb, Elisa - O14798

Item Number: CM0016059

Introduction

Rabbit polyclonal antibody targeting human TNFRSF10C (DcR1/TRAIL-R3). Suitable for Western blot and ELISA, with cross-reactivity to mouse and rat. Explores decoy receptor function in TRAIL-induced apoptosis.

[Learn More](#)

Parameter	Value
Product Name	TNFRSF10C Rabbit pAb
Remarks/Alias	LIT, DCR1, TRID, CD263, TRAILR3, TRAIL-R3, DCR1-TNFR, TNFRSF10C
Species	Human
GeneID (Human)	8794
GeneID	8794
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 150-238 of human TNFRSF10C (NP_003832.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O14798
Protein Weight	27kDa
Shipping	Ice bag