

Anti-Tnfrsf10D Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Ubn6

Item Number: CM0025368

Introduction

Rabbit polyclonal antibody against human TNFRSF10D (DcR2/TRAIL-R4). Validated for WB, IHC-P, ELISA. Detects endogenous TRAIL decoy receptor; cross-reacts with mouse. Ideal for apoptosis and cancer research.

[Learn More](#)

Parameter	Value
Product Name	TNFRSF10D Rabbit pAb
Remarks/Alias	DCR2; CD264; TRUND; TRAILR4; TRAIL-R4; TNFRSF10D
Species	Human
Gene ID (Human)	8793
Gene ID	8793
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 257-386 of human TNFRSF10D (NP_003831.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UBN6
Protein Weight	42kDa
Shipping	Ice bag