

Anti-Tnfrsf10D/Trailr4 Monoclonal Antibody For Fc - Q9Ubn6

Item Number: CM0007240

Introduction

This rabbit monoclonal antibody against human TNFRSF10D/TRAILR4 (DcR2, TRAIL-R4) is produced using a recombinant protein immunogen and validated for FC. Suitable for research on TRAIL-induced apoptosis and decoy receptor function.

[Learn More](#)

Parameter	Value
Product Name	TNFRSF10D/TRAILR4 Rabbit mAb
Remarks/Alias	DCR2; CD264; TRUNDD; TRAILR4; TRAIL-R4; TNFRSF10D/TRAILR4
Species	Human
Gene ID (Human)	8793
Gene ID	8793
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 56-211 of human TNFRSF10D/TRAILR4 (NP_003831.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9UBN6
Protein Weight	42kDa
Shipping	Ice bag