


Anti-Tnfrsf11A Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Y6Q6

Item Number: CM0018228



Introduction

Rabbit polyclonal antibody targeting TNFRSF11A (RANK/CD265) validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Suitable for osteoclastogenesis and immune signaling research.

[Learn More](#)

Parameter	Value
Product Name	TNFRSF11A Rabbit pAb
Remarks/Alias	FEO; OFE; ODFR; OST5; PDB2; RANK; CD265; OPTB7; TRANCER; LOH18CR1; TRANCE-R; TNFRSF11A
Species	Human
Gene ID (Human)	8792
Gene ID	8792
Immunogen	Recombinant protein/Recombinant fusion protein containing a sequence corresponding to amino acids 466-616 of human TNFRSF11A (NP_003830.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y6Q6
Protein Weight	66kDa
Shipping	Ice bag