


Anti-Tnfrsf11B Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - O00300

Item Number: CM0020380



Introduction

TNFRSF11B rabbit polyclonal antibody for Western blot, IF-paraffin, IHC-paraffin, and ELISA. Recognizes human, mouse, and rat TNFRSF11B. Immunogen: recombinant protein corresponding to amino acids 22-280 of human TNFRSF11B (NP_002537.3). Research use: osteoclastogenesis, bone remodeling, and apoptosis.

[Learn More](#)

Parameter	Value
Product Name	TNFRSF11B Rabbit pAb
Remarks/Alias	OPG; TR1; OCIF; PDB5; TNFRSF11B
Species	Human
Gene ID (Human)	4982
Gene ID	4982
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-280 of human TNFRSF11B (NP_002537.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00300
Protein Weight	46kDa
Shipping	Ice bag