


Anti-Tnfrsf11B Rabbit Polyclonal Antibody For Wb And Elisa - O00300

Item Number: CM0013401



Introduction

High-quality rabbit polyclonal antibody targeting human TNFRSF11B (Osteoprotegerin, OPG). Validated for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for studying osteoclast differentiation and bone homeostasis.

[Learn More](#)

Parameter	Value
Product Name	TNFRSF11B Rabbit pAb
Remarks/Alias	OPG; TR1; OCIF; PDB5; TNFRSF11B
Species	Human
GeneID	4982
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-280 of human TNFRSF11B (NP_002537.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00300
Protein Weight	46kDa
Shipping	Ice bag