


Anti-Tnfrsf13C Rabbit Polyclonal Antibody For Wb And Elisa - Q96Rj3

Item Number: CM0013868



Introduction

Rabbit polyclonal antibody against human TNFRSF13C (BAFF-R, CD268) validated for WB and ELISA, cross-reacts with mouse. Ideal for B-cell survival and immunodeficiency research.

[Learn More](#)

Parameter	Value
Product Name	TNFRSF13C Rabbit pAb
Remarks/Alias	BAFFR; CD268; CVID4; BAFF-R; BROMIX; prolixin; TNFRSF13C
Species	Human
GeneID (Human)	115650
GeneID	115650
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-80 of human TNFRSF13C (NP_443177.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96RJ3
Protein Weight	19kDa
Shipping	Ice bag