

Anti-Tnfrsf17 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q02223

Item Number: CM0017768



Introduction

High-quality rabbit polyclonal antibody targeting human TNFRSF17 (BCMA), validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for B-cell biology and humoral immunity research.

[Learn More](#)

Parameter	Value
Product Name	TNFRSF17 Rabbit pAb
Remarks/Alias	BCM; BCMA; CD269; TNFRSF13A; TNFRSF17
Species	Human
GeneID (human)	608
GeneID	608
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human TNFRSF17 (NP_001183.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02223
Protein Weight	20kDa
Shipping	Ice bag