

Anti-Tnfrsf25 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - Q93038

Item Number: CM0014468



Introduction

Rabbit polyclonal antibody targeting human TNFRSF25 (DR3, LARD). Suitable for Western blot, immunofluorescence, and ELISA. Cross-reacts with mouse and rat. Detects 45kDa protein. For apoptosis and lymphocyte homeostasis studies.

[Learn More](#)

Parameter	Value
Product Name	TNFRSF25 Rabbit pAb
Remarks/Alias	DR3, TR3, DDR3, LARD, APO-3, TRAMP, WSL-1, GEF720, WSL-LR, PLEKHG5, TNFRSF12, TNFRSF25
Species	Human
Gene ID (Human)	8718
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-190 of human TNFRSF25 (NP_683866.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q93038
Protein Weight	45kDa
Shipping	Ice bag