


Anti-Tnfsf10 Rabbit Polyclonal Antibody For Wb, Elisa - P50591

Item Number: CM0020668



Introduction

High-quality TNFSF10 rabbit polyclonal antibody validated for WB and ELISA. Detects human TNFSF10/TRAIL, a key apoptosis-inducing cytokine, suitable for research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	TNFSF10 Rabbit pAb
Remarks/Alias	TL2; APO2L; CD253; TANCRC; TRAIL; Apo-2L; TNLG6A; TNFSF10
Species	Human
GeneID (Human)	8743
GeneID	8743
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P50591
Protein Weight	33kDa
Shipping	Ice bag