


Anti-Tnfsf11 Rabbit Polyclonal Antibody For Ihc-P, Elisa - O14788

Item Number: CM0021571



Introduction

Rabbit polyclonal antibody to TNFSF11 (RANKL, OPGL) for IHC-P and ELISA applications; reacts with human samples. Ideal for studies on osteoclast differentiation, immune regulation, and bone resorption.

[Learn More](#)

Parameter	Value
Product Name	TNFSF11 Rabbit pAb
Remarks/Alias	ODF; OPGL; sOdf; CD254; OPTB2; RANKL; TNLG6B; TRANCE; hRANKL2; TNFSF11
Species	Human
GeneID (Human)	8600
GeneID	8600
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 165-317 of human TNFSF11 (NP_003692.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	O14788
Protein Weight	35kDa
Shipping	Ice bag