


Anti-Tnfsf13 (April) Rabbit Polyclonal Antibody For Wb, Elisa - 075888

Item Number: CM0016523



Introduction

Anti-TNFSF13 (APRIL) rabbit polyclonal antibody for WB and ELISA. Recognizes human TNFSF13, a cytokine involved in tumor cell growth and immune processes. Ideal for cancer and immunology research.

[Learn More](#)

Parameter	Value
Product Name	TNFSF13 Rabbit pAb
Remarks/Alias	APRIL; CD256; TALL2; ZTNF2; TALL-2; TNLG7B; TRDL-1; UNQ383/PRO715; TNFSF13
Species	Human
GeneID (human)	8741
GeneID	8741
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 105-247 of human TNFSF13 (NP_742085.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	075888
Protein Weight	27kDa
Shipping	Ice bag