


Anti-Tnfsf13B Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Y275

Item Number: CM0028572



Introduction

High-quality rabbit polyclonal antibody against human TNFSF13B (BAFF). Validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for B-cell immunology research.

[Learn More](#)

Parameter	Value
Product Name	TNFSF13B Rabbit pAb
Remarks/Alias	DTL; BAFF; BLYS; CD257; TALL1; THANK; ZTNF4; TALL-1; TNLG7A; TNFSF20; TNFSF13B
Species	Human
GeneID (Human)	10673
GeneID	10673
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 90-160 of human TNFSF13B (NP_006564.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y275
Protein Weight	31kDa
Shipping	Ice bag