

Anti-Tnk1 Polyclonal Antibody For Wb And Elisa - Q13470

Item Number: CM0022054



Introduction

Rabbit polyclonal antibody against human non-receptor tyrosine kinase TNK1, suitable for WB and ELISA. Involved in tumor suppression and Ras-MAPK pathway regulation. Ideal for cancer and signal transduction research.

[Learn More](#)

Parameter	Value
Product Name	TNK1 Rabbit pAb
Remarks/Alias	KOS1; TNK1
Species	Human
GeneID (Human)	8711
GeneID	8711
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-380 of human TNK1 (NP_001238831.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q13470
Protein Weight	72 kDa
Shipping	Ice bag