

Anti-Tnks1Bp1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9C0C2

Item Number: CM0022735

Introduction

Rabbit polyclonal antibody targeting human TNKS1BP1 (TAB182), validated for WB, IHC-P, and ELISA, with cross-reactivity to mouse and rat. Suitable for chromatin, cytoskeleton, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	TNKS1BP1 Rabbit pAb
Remarks/Alias	TAB182; TNKS1BP1
Species	Human
GeneID (Human)	85456
GeneID	85456
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 700-1000 of human TNKS1BP1 (NP_203754.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9C0C2
Protein Weight	182kDa
Shipping	Ice bag