


Anti-Tnnc1 Rabbit Monoclonal Antibody For Wb And Elisa - P63316

Item Number: CM0005379



Introduction

TNNC1 Rabbit mAb is a monoclonal antibody against troponin C, validated for Western blot and ELISA, with cross-reactivity to mouse and rat. Ideal for cardiac muscle research and disease studies.

[Learn More](#)

Parameter	Value
Product Name	TNNC1 Rabbit mAb
Remarks/Alias	TNC; TN-C; TNNC; CMD1Z; CMH13; TNNC1
Species	Human
Gene ID	7134 (Human)
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human TNNC1 (P63316)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P63316
Protein Weight	18kDa
Shipping	Ice bag