


Anti-Tnnc1 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P63316

Item Number: CM0022960



Introduction

Rabbit polyclonal antibody against human TNNC1 (Troponin C, slow skeletal and cardiac muscles). Validated for WB, IF-P, ELISA. Cross-reacts with mouse and rat. Suitable for muscle contraction and cardiomyopathy research.

[Learn More](#)

Parameter	Value
Product Name	TNNC1 Rabbit pAb
Remarks/Alias	TNC; TN-C; TNNC; CMD1Z; CMH13; TNNC1
Species	Human
GeneID (Human)	7134
GeneID	7134
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human TNNC1 (NP_003271.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P63316
Protein Weight	18kDa
Shipping	Ice bag