

Anti-Tnnc2 Polyclonal Antibody For Wb, Ihc-P, Elisa - P02585

Item Number: CM0028864



Introduction

TNNC2 rabbit polyclonal antibody validated for Western blot (WB), immunohistochemistry (IHC-P), and ELISA in human, mouse, and rat samples. Recombinant immunogen spanning amino acids 1-160 of human troponin C, skeletal muscle. Suitable for muscle contraction research.

[Learn More](#)

Parameter	Value
Product Name	TNNC2 Rabbit pAb
Remarks/Alias	FAP85; CFAP85; CMYP15; MYONRI; TNNC2
Species	Human
GeneID (Human)	7125
GeneID	7125
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human TNNC2 (NP_003270.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02585
Protein Weight	18kDa
Shipping	Ice bag