

Anti-Tnni3K Polyclonal Antibody For Wb, Elisa - Q59H18

Item Number: CM0030422



Introduction

High-quality rabbit polyclonal antibody against human TNNI3K. Validated for Western blot and ELISA. Cross-reacts with human, mouse, and rat. Ideal for cardiac function and cardiomyopathy research.

[Learn More](#)

Parameter	Description
Product Name	TNNI3K Rabbit pAb
Remarks/Alias	CARK; CCDD; TNNI3K
Species	Human
GeneID (Human)	51086
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-75 of human TNNI3K (NP_057062.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q59H18
Protein Weight	93kDa
Shipping	Ice bag