

Anti-Tnnt1 Polyclonal Antibody For Wb, Elisa - P13805

Item Number: CM0014150

Introduction

Rabbit polyclonal antibody against human TNNT1 (slow skeletal troponin T), validated for Western blot and ELISA. Suitable for muscle contraction and nemaline myopathy research.

[Learn More](#)

Parameter	Value
Product Name	TNNT1 Rabbit pAb
Remarks/Alias	ANM; TNT; NEM5; STNT; TNTS; TNNT1
Species	Human
GeneID	7138
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-262 of human TNNT1 (NP_001119604.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P13805
Protein Weight	33kDa
Shipping	Ice bag