

## Anti-Tnpo1 Polyclonal Antibody For Wb, Elisa - Q92973

Item Number: CM0025842



### Introduction

Rabbit polyclonal antibody against human TNPO1 (Transportin-1), validated for Western blot and ELISA. Shows cross-reactivity with human, mouse, and rat. Ideal for nuclear transport and importin studies.

[Learn More](#)

| Parameter        | Value                                                                                                 |
|------------------|-------------------------------------------------------------------------------------------------------|
| Product Name     | TNPO1 Rabbit pAb                                                                                      |
| Remarks/Alias    | MIP; TRN; IPO2; MIP1; KPNB2; TNPO1                                                                    |
| Species          | Human                                                                                                 |
| Gene ID (Human)  | 3842                                                                                                  |
| Gene ID          | 3842                                                                                                  |
| Immunogen        | Synthetic peptide corresponding to a sequence within amino acids 790-890 of human TNPO1 (NP_002261.3) |
| Source           | Rabbit                                                                                                |
| Category         | Polyclonal Antibodies                                                                                 |
| Application      | WB, ELISA                                                                                             |
| Cross Reactivity | Human, Mouse, Rat                                                                                     |
| SWISS            | Q92973                                                                                                |
| Protein Weight   | 102kDa                                                                                                |
| Shipping         | Ice bag                                                                                               |