

Anti-Tnr Polyclonal Antibody For Wb, Elisa - Q92752

Item Number: CM0018152



Introduction

Rabbit polyclonal antibody against human tenascin-R (TNR), validated for WB and ELISA with cross-reactivity to mouse and rat. Suitable for neural extracellular matrix research.

[Learn More](#)

Parameter	Value
Product Name	TNR Rabbit pAb
Remarks/Alias	TN-R; NEDSTO; TNR
Species	Human
GeneID (Human)	7143
GeneID	7143
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1059-1358 of human TNR (NP_003276.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92752
Protein Weight	150kDa
Shipping	Ice bag