

Anti-Tnrc6A Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Ndv7

Item Number: CM0026390



Introduction

Anti-TNRC6A rabbit polyclonal antibody, suitable for WB, IF/ICC, and ELISA, cross-reacts with human, mouse, and rat. Ideal for studying miRNA-mediated gene silencing and P-body components.

[Learn More](#)

Parameter	Value
Product Name	TNRC6A Rabbit pAb
Remarks/Alias	GW1; FAME6; GW182; TNRC6; CAGH26; TNRC6A
Species	Human
Gene ID (Human)	27327
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NDV7
Protein Weight	210kDa
Shipping	Ice bag