

Anti-Tob2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q14106

Item Number: CM0028275

Introduction

Rabbit polyclonal antibody targeting TOB2, validated for WB, IHC-P, ELISA with cross-reactivity to human, mouse, and rat. TOB2 is an anti-proliferative protein that inhibits G0/G1 to S phase transition, localized in cytoplasm. Ideal for cell cycle and oncology research.

[Learn More](#)

Parameter	Value
Product Name	TOB2 Rabbit pAb
Remarks/Alias	TOB4, TOBL, APRO5, TROB2, TOB2
Species	Human
Gene ID (Human)	10766
Gene ID	10766
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-70 of human TOB2 (NP_057356.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14106
Protein Weight	37kDa
Shipping	Ice bag