


Anti-Tollip Polyclonal Antibody For Wb, Elisa - Q9H0E2

Item Number: CM0022254



Introduction

High-quality rabbit polyclonal antibody targeting human TOLLIP, validated for Western blot and ELISA. Suitable for human, mouse, rat samples; ideal for innate immunity and autophagy research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] TOLLIP Rabbit pAb
Remarks/Alias	IL-1RAcPIP; [KD Validated] TOLLIP
Species	Human
GeneID (Human)	54472
GeneID	54472
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H0E2
Protein Weight	30kDa
Shipping	Ice bag