


Anti-Tollip Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9H0E2

Item Number: CM0021395



Introduction

TOLLIP Rabbit pAb is a polyclonal antibody sourced from rabbit, validated for WB, IHC-P, IF/ICC, ELISA, and reacts with human, mouse, rat. It targets Toll-interacting protein, a key regulator of IL-1 and TLR signaling and autophagy.

[Learn More](#)

Parameter	Value
Product Name	TOLLIP Rabbit pAb
Remarks/Alias	IL-1RAcPI; TOLLIP
Species	Human
GeneID (Human)	54472
GeneID	54472
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-274 of human TOLLIP (NP_061882.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H0E2
Protein Weight	30kDa
Shipping	Ice bag