


Anti-Tollip Rabbit Monoclonal Antibody (Kd Validated) For Wb, Ihc-P, Elisa - Q9H0E2

Item Number: CM0007699



Introduction

High-quality rabbit monoclonal antibody against TOLLIP, KD validated, for Western blot, IHC-P, and ELISA. Cross-reacts with human, mouse, rat. Useful for studying innate immunity, Toll-like receptor signaling, and autophagy pathways.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] TOLLIP Rabbit mAb
Remarks/Alias	TOLLIP; IL-1RAcP; toll-interacting protein
Species	Human
Gene ID	54472
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-274 of human TOLLIP (NP_061882.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H0E2
Protein Weight	22kDa/30kDa
Shipping	Ice bag