


Anti-Tom1 Rabbit Polyclonal Antibody For Wb And Elisa - O60784

Item Number: CM0026924



Introduction

Rabbit polyclonal antibody against human TOM1, validated for WB and ELISA. Cross-reacts with mouse and rat. Target of Myb1 involved in endosomal trafficking of ubiquitinated proteins. Ideal for autophagy and signaling studies.

[Learn More](#)

Parameter	Value
Product Name	TOM1 Rabbit pAb
Remarks/Alias	IMD85; TOM1
Species	Human
Gene ID (Human)	10043
Gene ID	10043
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 314-493 of human TOM1 (NP_001129204.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O60784
Protein Weight	54kDa
Shipping	Ice bag