


Anti-Tom1L1 Rabbit Polyclonal Antibody For Wb, Elisa - 075674

Item Number: CM0022076



Introduction

Anti-TOM1L1 Rabbit Polyclonal Antibody validated for Western blot and ELISA, reacting with human, mouse, and rat TOM1L1. TOM1L1 is an adapter protein involved in signaling pathways, localized to Golgi and endosomes.

[Learn More](#)

Parameter	Value
Product Name	TOM1L1 Rabbit pAb
Remarks/Alias	SRCASM; TOM1L1
Species	Human
GeneID (human)	10040
GeneID	10040
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 250-350 of human TOM1L1 (NP_005477.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075674
Protein Weight	53kDa
Shipping	Ice bag