

Anti-Tom1L2 Polyclonal Antibody For Wb, Elisa - Q6Zvm7

Item Number: CM0024515

Introduction

TOM1L2 Rabbit pAb is a polyclonal antibody for WB and ELISA, detecting human, mouse, and rat TOM1L2 (~56 kDa). Useful for research on endocytic trafficking and mitogenic signaling.

[Learn More](#)

Parameter	Value
Product Name	TOM1L2 Rabbit pAb
Remarks/Alias	TOM1L2
Species	Human
GeneID (Human)	146691
GeneID	146691
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 310-450 of human TOM1L2 (NP_001076437.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6ZVM7
Protein Weight	56kDa
Shipping	Ice bag