

Anti-Tom20 Monoclonal Antibody (Conjugated) For Wb And If - Q15388

Item Number: CM0010422



Introduction

Pre-conjugated Anti-TOM20 rabbit monoclonal antibody for direct detection of human TOM20 in WB, IF/ICC, and IF-P. Validated in human, mouse, and rat samples. Ideal for mitochondrial import and mitophagy research.

[Learn More](#)

Parameter	Value
Product Name	Conjugated TOM20 Rabbit mAb
Remarks/Alias	MAS20; MOM19; TOM20
Species	Human
Gene ID	9804
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-145 of human TOM20 (NP_055580.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15388
Protein Weight	16kDa
Shipping	Ice bag