

Anti-Tom20 Monoclonal Antibody For Wb, Ip, If/Icc, If-P, Ihc-P, Elisa - Q15388

Item Number: CM0005256



Introduction

Mouse monoclonal antibody against human TOM20, validated for WB, IP, IF/ICC, IF-P, IHC-P, and ELISA. Detects human, mouse, rat. Ideal for mitochondrial import studies.

[Learn More](#)

Parameter	Specification
Product Name	TOM20 Mouse mAb
Remarks/Alias	MAS20; MOM19; TOM20
Species	Human
GeneID (human)	9804
GeneID	9804
Immunogen	Recombinant protein of human TOM20.
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15388
Protein Weight	16 kDa
Shipping	Ice bag