

Anti-Tom20 Polyclonal Antibody For Wb, Ip, If/Icc, Ihc-P, Elisa - Q15388

Item Number: CM0030288



Introduction

TOM20 Rabbit pAb is a polyclonal antibody recognizing human, mouse, and rat TOM20. Suitable for WB, IP, IF/ICC, IHC-P, ELISA, it is ideal for mitochondrial import, mitophagy, and dynamics studies.

[Learn More](#)

Parameter	Value
Product Name	TOM20 Rabbit pAb
Remarks/Alias	MAS20; MOM19; TOM20
Species	Human
Gene ID (Human)	9804
Gene ID	9804
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-145 of human TOM20 (NP_055580.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15388
Protein Weight	16 kDa
Shipping	Ice bag