

Anti-Tom20 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q15388

Item Number: CM0002726



Introduction

High-quality rabbit monoclonal antibody targeting human TOM20, validated for WB, IHC-P, IF/ICC, IP, and ELISA. Detects human, mouse, and rat TOM20, a key mitochondrial import receptor. Ideal for mitochondrial research and mitophagy studies.

[Learn More](#)

Parameter	Value
Product Name	TOM20 Rabbit mAb
Remarks/Alias	MAS20; MOM19; TOM20
Species	Human
GeneID (human)	9804
GeneID	9804
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-145 of human TOM20 (NP_055580.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15388
Protein Weight	16kDa
Shipping	Ice bag