

Anti-Tom20 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q15388

Item Number: CM0019310



Introduction

TOM20 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA in Human, Mouse, Rat. Targets mitochondrial import receptor TOM20 homolog (Q15388) functioning in mitophagy and protein translocation. Ideal for mitochondrial research.

[Learn More](#)

Parameter	Value
Product Name	TOM20 Rabbit pAb
Remarks/Alias	MAS20; MOM19; TOM20; 20
Species	Human
GeneID	9804
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-145 of human TOM20 (NP_055580.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15388
Protein Weight	16kDa
Shipping	Ice bag