

Anti-Tom20 Rat Monoclonal Antibody For Wb, If, Ihc, Elisa - Q15388

Item Number: CM0007800



Introduction

High-specificity anti-TOM20 rat monoclonal antibody validated for WB, IF/ICC, IF-P, IHC-P, and ELISA. Detects human, mouse, and rat TOM20, a mitochondrial import receptor. Ideal for mitochondrial function and mitophagy research.

[Learn More](#)

Parameter	Value
Product Name	TOM20 Rat mAb
Remarks/Alias	MAS20; MOM19; TOM20
Species	Human
Gene ID (Human)	9804
Gene ID	9804
Immunogen	Recombinant protein of human TOMM20.
Source	Rat
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15388
Protein Weight	16kDa
Shipping	Ice bag